

# **Ballistic Pendulum**

## Purpose

Use the ballistic pendulum to observe the collision phenomenon of two objects and record their results. The experiment can examine the conservation of momentum and kinetic energy. In normal situation, a collision can be classified into two types, perfectly elastic collision and perfectly inelastic collision. In this experiment, we will discuss the situation of perfectly inelastic collision.

### Theory

From the conservation of momentum, we know that

$$\Delta \vec{p} = m \Delta \vec{v} = \vec{F} \Delta t$$

When two objects oscillate, the motion system conserves the momentum. Based on the law of conservation of energy, the system's energy can be classified into two types:

1. The energy is completely transferred to kinetic energy which is conservation of kinetic energy.

2. If the energy cannot be completely transferred to kinetic energy due to fraction, resistance or deformation, it loses its kinetic energy or is transferred to other form of energy such as heat but the total energy in the system is conserved. The kinetic energy is not conserved in the collision.

Suppose masses of two objects are  $m_1$  and  $m_2$ , original velocities of objects are  $v_1$  and  $v_2$ , and final velocities of objects are  $u_1$  and  $u_2$ . Based on the law of conservation of momentum, the formula can be expressed as:

$$m_1 \, \vec{v}_1 + m_2 \, \vec{v}_2 = m_1 \, \vec{u}_1 + m_2 \, \vec{u}_2$$

The velocity difference before the collision of two objects is  $v_1 - v_2$ . The velocity difference after collision is  $u_1 - u_2$ . The ratio of the velocity differences is the coefficient of restitution e.



$$e = \frac{u_2 - u_1}{v_1 - v_2}$$

According to above description, if we ignore the ground fraction and the air resistance in a motion system, the collision can be classified into

1. Perfectly elastic collision (e = 1)









2. Perfectly inelastic collision (e = 0)



3. Non-completely elastic collision (0 < e < 1)

The ballistic pendulum in this experiment is perfectly inelastic collision.



### **Ballistic Pendulum**



From the law of conservation of energy and conservation of momentum, we know that

 $P_i = P_f$  and mv = (M + m) V(K.E)  $_{bottom} = (P.E)_{top}$  $\frac{1}{2}(M + m)V^2 = (M + m)gh$  $h = L \times (1 - \cos \theta)$ 

Total energy may or may not be lost during the experiment process. In this experiment, we assume the lost of energy is small and can be ignored.



### Instruments

| Instrument List |   |     |     |  |     |
|-----------------|---|-----|-----|--|-----|
| No.             | Accessory   | Qty | No. | Accessory  | Qty |
| 1.              | Multi-function aluminum base  | 1   | 2.  | Launcher holder (with                                    | 1   |
|                 | (including slide implement *1)  |     |     | large knob* 4)   |     |
| 3.              | Projectile launcher   | 1   | 4.  | Set of steel balls                                       | 1   |
|                 | (including slide implement*1)   |     |     | (D16mm×3)  |     |
| 5.              | Loading rod   | 1   | 6.  | Level  | 1   |
| 7.              | Accessory of ballistic<br>pendulum (including ballistic<br>pendulum, angle disk,<br>indicator, fixed slotting and<br>board) | 1   | 8.  | Movable connector base                                   | 1   |
| 9.              | Weight 50g  | 2   | 10. | Set of photogate fixer<br>(with clips)                   | 1   |
| 11.             | Photogate electronic timer<br>(with power supply 12 VDC)<br>Note: Additional purchase                                       |     | 12. | Photogate (with metal<br>fixed knob)<br>Note: Additional |     |
|                 | Note: Additional purchase   |     |     | Note: Additional purchase                                |     |







Atis Scientific Instruments Co.,Ltd Address: 1F., No.18, Nanming St., South Dist., Tainan City 702, Taiwan (R.O.C.) E-mail:atis@atissi.com.tw Tel: (886) -6-2925201 Fax: (886) -6-2611476 Mobile:+886-9-8006-1128 Website: www.atis.com.tw

All rights reserved